



COVIDIEN ENERGY.

Performance Amplified.



COVIDIEN

LigaSure™ Maryland Jaw Open/Laparoscopic Sealer/Divider

Reach for perfection.

Efficient. Versatile. Multifunctional.



LigaSure™ Maryland jaw combines one-step sealing with the functionality of a Maryland dissector, atraumatic grasper and cold scissors with the reliability of LigaSure™ technology.

EFFICIENT

- One-step sealing provides efficient transection speed¹
- Design allows for minimal steps when sealing and dividing¹
- The actions to grasp, seal and cut are simple and intuitive¹

VERSATILE

- Enhanced blunt dissection¹
- Improved tip visualization¹
- Reduced instrument exchanges¹

MULTIFUNCTIONAL

- A Maryland dissector, grasper and cold scissors all from one device
- Consistency, control and safety of LigaSure™ technology

One-step sealing.



GRASPING AND MANIPULATING TISSUE

Place the tissue in the jaws and pull back on the handle. The **first click** indicates the end of the grasp zone and alerts the user that additional pressure will activate energy.



SEALING VESSELS AND TISSUE BUNDLES

Deliver RF energy by squeezing the lever until the **purple activation button clicks**. Continue to hold the lever closed until the seal cycle is complete.*

*A continuous tone sounds to indicate that the vessel or tissue is being sealed. When the seal cycle is complete, a short end tone sounds (two consecutive tones) and RF output ceases.



DIVIDING TISSUE

Pull the **cutting trigger**. Release the cutting trigger to retract the blade.



LigaSure™ technology is supported by over 300 peer-reviewed clinical studies and over eight million procedures.

WHAT SURGEONS HAVE SAID ABOUT THE LIGASURE™ MARYLAND DEVICE

After using the LigaSure™ Maryland device:

- 100% of general surgeons agreed that it will reduce instrument exchanges in surgery¹
- 90% of general surgeons agreed that the integrated cutter will reduce the need for additional cold scissors¹
- 80% of general surgeons agreed that the device will reduce the need for additional dissecting instruments¹
- The majority of surgeons agreed that the secure and atraumatic grasping will reduce the need for an additional grasper¹
- 100% of surgeons surveyed after using the instrument believe that the LigaSure™ Maryland device provides efficiency throughout the procedure¹

WHY MORE SURGEONS CHOOSE LIGASURE™ TECHNOLOGY

Compared to mechanical ligation techniques, LigaSure™ technology has been shown to significantly reduce:

- Operative blood loss in colorectal, gynecologic and urologic surgery⁴⁻⁹
- Perioperative blood transfusions in gynecologic, urologic and general surgery^{8,10,11}
- Procedure time in colorectal, gynecologic and urologic surgery^{4,6,7,9,12,13}
- Length of hospital stay in gynecologic and urologic surgery^{6,13}

Compared to other energy-based modalities, LigaSure™ technology has been shown to significantly reduce:

- Operative blood loss in colorectal and gynecologic surgery¹⁴⁻¹⁷
- Procedure time in colorectal and gynecologic surgery^{11, 17}



COVIDIEN LigaSure Maryland 5mm-37cm

Multiple lengths provide standardization across procedures²

- 23 cm provides access in deep and confined spaces¹
- 44 cm provides additional reach²

WIDE JAW APERTURE
13 mm

JAW LENGTH
20 mm

SHAFT LENGTH
23 cm, 37 cm, 44 cm

SHAFT ROTATION
350 Degrees

CUT LENGTH
18 mm

ACTIVATION METHOD
Purple Button

ENERGIZED BY
ForceTriad™ energy platform
(software version 3.60 and higher)

Ordering Information

Catalog Number	Quantity Per Package	Product Description
LF1723	6 each	LigaSure™ Maryland Device – 23 cm
LF1737	6 each	LigaSure™ Maryland Device – 37 cm
LF1744	6 each	LigaSure™ Maryland Device – 44 cm

1. Based on independent surgeon feedback collected during Covidien-sponsored labs conducted April 16-18, 2013 and April 30-May 3, 2013
 2. LigaSure™ Maryland device IFU
 3. Compared to standard length laparoscopic devices

REFERENCES

1. Based on independent surgeon feedback collected during Covidien-sponsored labs conducted April 16-18, 2013 and April 30-May 3, 2013.
2. LigaSure™ Maryland device IFU.
3. Compared to standard length laparoscopic devices.
4. Targarona, et al. A prospective randomized comparison of conventional electro-surgery, bipolar computer-controlled electro-surgery and ultrasonic dissection. Operative Outcome and Cost analysis. Surgical Innovation. 2005 Dec; 12(4): 339-344.
5. Manouras, et al. Sutureless open low anterior resection with total mesorectal excision for rectal cancer with the use of the electrothermal bipolar vessel sealing system. 2007. Med Sci Monit; 13(5): CR224-230.
6. Ding Z, et al. Use of LigaSure™ bipolar diathermy system in vaginal hysterectomy. Journal of Obstetrics Gynaecology. 2005 Jan; 25(1): 49-51.
7. Levy B, Emery L. Randomized Trial of Suture Versus Electrosurgical Bipolar Vessel Sealing in Vaginal Hysterectomy. Obstetrics and Gynecology. 2003 Jul; 102(1): 147-151.
8. Tamussino K, et al. Electrosurgical bipolar vessel sealing for radical abdominal hysterectomy. Gynecologic Oncology. 2005 Feb; 96(2): 320-322.
9. Leonardo C, et al. Laparoscopic nephrectomy using LigaSure™ system: preliminary experience. Journal of Endourology. 2005 Oct; 19(8): 976-8.
10. Daskalopoulos G et al. Electrothermal bipolar coagulation for radical prostatectomies and cystectomies: a preliminary case-controlled study. International Urology and Nephrology. 2004; 36(2): 181-185.
11. Romano F, et al. Laparoscopic Splenectomy: LigaSure™ versus EndoGIA: A Comparative Study. J Laparoendoscopic & Adv Surg Techniques. 2007. Vol. 17, No.6.
12. Cronje HS, et al. Electrosurgical bipolar vessel sealing during vaginal hysterectomy. Int J Gynaecol Obstet. 2005 Dec; 91(3): 243-5.
13. Metzelder ML, et al. Laparoscopic nephroureterectomy in children: a prospective study on LigaSure™ versus Clip/Ligation. Eur J Pediatr Surg. 2006 Aug; 16(4): 241-4.
14. Araki Y, et al. Clipless hand-assisted laparoscopic total colectomy using LigaSure™ Atlas. Kurume Medical Journal. 2004; 51(2): 105-8.
15. Campagnacci R, et al. Electrothermal bipolar vessel sealing device vs. ultrasonic coagulating shears in laparoscopic colectomies: a comparative study. Surg Endosc. 2007 Feb 8.
16. Takada et al. Comparative study of electrothermal bipolar vessel sealer and ultrasonic coagulating shears in laparoscopic colectomy. Surgical Endoscopy. 2005; 19:226-228.
17. Demirturk F, et al. Comparison of the use of electrothermal bipolar vessel sealer with harmonic scalpel in total laparoscopic hysterectomy. J Obst Gynaecol Res. 2007; 33 (3):341-345.



Change the energy in your OR. Contact your Covidien Sales Representative at 1-(800)255-8522 or visit Covidien.com/MarylandJaw.