



COVIDIEN

positive results for life™

WarmTouch™ Convective Air Warming System Not Too Hot. Not Too Cold.

WarmTouch WARMING BLANKETS

Covidien's WarmTouch warming blankets feature multi-layer, quilted construction. Strong two-ply material resists fluids, tears and punctures. Plus, a soft fabric outer layer protects patients from direct heat. Our special quilting prevents billowing, helps the blanket conform to the patient's body and promotes airflow. Quilting is also why WarmTouch blankets are more comfortable than tube-type blankets. They come in the following styles:

- Upper Body
- Lower Body
- Full Body
- Pediatric (2 sizes)
- Sterile Cardiac



THIS ONE IS JUST RIGHT

Cold compromises the immune systems of surgery patients. The WarmTouch system from Covidien keeps patients warm, which minimizes post-surgical complications and recovery times.

Our comfortable, durable blankets and our quiet, efficient warming unit give clinicians a convenient, cost-effective way to maintain patient temperatures—before, during and after surgery.

WARMTOUCH 5300A WARMING UNIT

Because thermal stability promotes patient stability, the 5300A warming unit warms up quickly (to 38°C in 30 seconds). It also has a boost setting for patients who need fast, short-term warming.

The 5300A has the highest airflow (55 CFM) of any convection warming system. It also features a HEPA filter to reduce airborne particulates. For further convenience, the 5300A has multiple mounting options—cart, bedrail, IV pole and floor—and a digital hour meter.

Find out today how the warmth and comfort of the WarmTouch system can help improve your surgical outcomes. Call 800-635-5267 or log on to www.covidien.com.

SPECIFICATIONS

Adult Warming Blankets

Size Uninflated: **Upper Body:** 28 in W x 82 in L
(71 cm W x 208 cm L)
Lower Body: 41 in W x 56 in L
(104 cm W x 142 cm L)
Full Body: 40 in W x 79 in L
(102 cm W x 165 cm L)

Size Inflated: **Upper Body:** 24 in W x 72 in L
(61 cm W x 183 cm L)
Lower Body: 34 in W x 52 in L
(86 cm W x 132 cm L)
Full Body: 33 in W x 74 in L
(84 cm W x 188 cm L)

Pediatric Warming Blankets

Small Pediatric: 25 in W x 41 in L
(63 cm W x 104 cm L)
Large Pediatric: 30 in W x 57 in L
(76 cm W x 145 cm L)

Small Pediatric: 22 in W x 35 in L
(56 cm W x 89 cm L)
Large Pediatric: 27 in W x 49 in L
(69 cm W x 124 cm L)

Sterile Cardiac Warming Blanket

Cardiac: 40 in W x 65 in L
(102 cm W x 165 cm L)

Cardiac: 34 in W x 53 in L
(86 cm W x 135 cm L)

Construction: 2-ply material consisting of polyethylene film inner layer and nonwoven outer layers (latex free).

WARMTOUCH 5300A WARMING UNIT

Size: 15 in H x 16 in W x 11 in D
(38 cm x 41 cm x 28 cm)
Weight: 15 lb (7.3 kg)
Hose: 5 ft (1.5 m); collapsible and self-supporting
AC Power
Cord: 14-ft long (4.3 m); hospital grade
CSA Listing: CSA 22.2 No. 601.1-M90
UL Standard: UL544, Hospital Medical Equipment

AC Power: 115 VAC, 60 Hz
Current Draw: 10 amperes, maximum
Heater Size: 1000 watts
Warm-Up Time: Air temperature at nozzle exceeds
38°C in approximately 30 seconds
Air Flow Rate: Approximately 55 CFM
Noise Level: Less than 55 db
(normal conversation level)

ORDERING INFORMATION

DESCRIPTION	QUANTITY	CATALOG NUMBER
Upper Body Blanket with Head Drape	12/case	5030870
Lower Body Blanket	12/case	5030880
Pediatric Blanket	12/case	5030850
Cardiac Blanket (Sterile)	12/case	5030860
PACU Blanket, Adult	25/case	5030810
PACU Blanket, Pediatric	12/case	5030840
WarmTouch 5300A Warming Unit	1/box	5015300A
Replacement Filter for WarmTouch 5300A Warming Unit	1/box	502-2200

For more information, please contact your local sales representative or visit our website at www.covidien.com



COVIDIEN, COVIDIEN with Logo and TM marked brands are trademarks of Covidien AG or an affiliate.
© 2008 Covidien AG or an affiliate. All rights reserved.
AP.a 3982v1 WT06507