

# Dakota™

Nitinol Stone Retrieval Device  
with OpenSure™ Handle

**EVERY DETAIL MATTERS**

When it comes to removing complex stones of varying size and composition in challenging locations,


# EVERY DETAIL MATTERS.

## LAST LONGER.



Capture, reposition, and remove multiple stones and fragments, all in a single procedure.

### Durability of Stone Retrieval Devices

Boston Scientific's Dakota™ Retrieval Device versus Cook's NGage™ Stone Extractor



In a Bench Study, Dakota's Retrieval Device demonstrated to be more durable than the NGage Stone Extractor when capturing larger stones ranging from 8mm – 10mm. **100% of the Dakota Retrieval Device survived 75 cycles without breakage compared to 23% of the NGage retrieval device.**<sup>1</sup>





## RELEASE LARGER.

Quickly release large, complex kidney stones on demand with a single press of the patented OpenSure™ Handle:


- **8mm device opens 50% larger**
- **11mm device opens 39% larger<sup>1</sup>**

## REACH DEEPER.

Our flexible, three-arm basket opens and closes at full deflection to reach the lower pole of the kidney in patients of tortuous anatomy.

## CAPTURE SMALLER.

A high-performance nitinol hybrid grasper enables head-on capture of small 1mm stones.



## Ordering information

Product Code	Description	Quantity
M0063905000	Dakota 1.9F x 8mm x 120cm	Single Unit
M0063905010	Dakota 1.9F x 11mm x 120cm	Single Unit

Every detail of this break-through hybrid basket is designed to bring you one step closer to your ultimate goal – **increasing the number of stone-free patients and improving outcomes.**

**WE'VE LEFT NO STONE UNTURNED.**

**Please consult your sales representative for more information and ordering details.**


**To learn more visit  
[www.bostonscientific.com/Dakota](http://www.bostonscientific.com/Dakota)**

1 Bench Test results may not necessarily be indicative of clinical performance. Data on file.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.



**Boston  
Scientific**  
Advancing science for life™

Boston Scientific Corporation  
300 Boston Scientific Way  
Marlborough, MA 01752  
[www.bostonscientific.com](http://www.bostonscientific.com)

© 2016 by Boston Scientific Corporation  
or its affiliates. All rights reserved.

URO-434806-AA OCT 2016